



# CONCEPT LIFE SCIENCES

DELIVERING SCIENCE

## DATA SHEET

## CYP 450 INHIBITION

Inhibitory potency ( $IC_{50}$ ) provides invaluable information regarding drug interaction potential. The  $IC_{50}$  of test compound against the five most important drug metabolising enzymes is determined using specific drug probe substrates. Analysis is by LC-MSMS.

**Deliverable:**  $IC_{50}$  ( $\mu M$ ) of test compound against CYP1A2, CYP2C9, CYP2C19, CYP2D6 & CYP3A4.

### CUSTOMER PROVIDES

Compound identifier and molecular formula.

Test: 70 $\mu$ L of 10mM in DMSO or 0.5mg solid.

### ENZYMES

Pooled recombinant E.coli membranes expressing human CYPs.

Table 1. Substrates & standard inhibitors:

CYP	SUBSTRATE AT $K_m$	METABOLITE PRODUCED	STANDARD INHIBITOR
1A2	Phenacetin	Acetaminophen	$\alpha$ -Naphthoflavone
2C9	Piclofenac	4'OH Diclofenac	Sulfaphenazole
2C19	S-mephenytoin	4'OH Mephenytoin	Tranylcypromine
2D6	Bufuralol	1'OH Bufuralol	Quinidine
3A4	Midazolam	1'OH Midazolam	Ketoconazole

### TEST COMPOUND

Concentration range 0, 0.15, 0.5, 1.5, 5, 15 & 50 $\mu$ M.

### FORMAT

96-well plate, 200 $\mu$ L incubation volume.

### PROTOCOL

CYPs are pooled in a ratio such that biotransformation of each probe substrate is specific for the particular CYP450. Substrates (at their  $K_m$ ) are also pooled within the same solution to create an enzyme-substrate stock, which is aliquoted into to each well of the microplate. Test compound is added to appropriate wells. The final solvent concentration is 1.0% DMSO. After equilibration to 37°C, addition of NADPH and mixing initiates substrate biotransformation. Each concentration of test compound is assayed against five CYPs in the same well simultaneously. The incubation is stopped at  $t = 10$ min by removal of the plate from the shaking incubator, followed by the addition of 200 $\mu$ L of ice cold methanol containing internal standard, and mixing.

The quenched samples are stored at -20°C for a minimum of 4hr, then centrifuged (3800rpm at 4°C for 20min) to precipitate the protein.

### QUANTITATION

The supernatants analysed by LC-MS/MS using Concept Life Sciences generic analytical methods to measure metabolite formation.





# CONCEPT LIFE SCIENCES

DELIVERING SCIENCE

## DATA ANALYSIS AND RESULTS

The enzyme activity within each test compound incubation is compared with Control (contains solvent but no test compound). For each concentration of test:

$$\% \text{ Control activity} = \frac{\text{Test Peak Area}}{\text{Internal Standard Peak Area in Test}} \times \frac{\text{Internal Standard Peak Area in Control}}{\text{Control Peak Area}} \times 100\%$$

% control activities obtained from the above calculation are plotted against test concentration.

IC<sub>50</sub> = test compound concentration producing 50% control activity (=50% inhibition).


## +VE CONTROLS

Control compound, one for each CYP, is used in every assay run (see bar chart below for identities).


## RESULTS

6-point inhibition curves are used to calculate IC<sub>50</sub>'s for all five CYPs.

### 3A4 KETOCONAZOLE INHIBITION



### 2D6 QUINIDINE INHIBITION



## CONCEPT IC50

